

The **BABY LEO VITA 48** allows all new borns to get the ideal temperature and warms up the « low weight » ones non able to maintain their necessary normothermia the first days of their life. BABY LEO VITA 48 warming cot is dedicated to neonatology and maternity departments.



- Easy to use
- Comfortable and safe for the baby
- Made in France !



BABY LEO VITA 48 Warming Cot

The **BABY LEO VITA 48** is a compromise between the incubator and the ordinary cot. This allows a link between those 2 universes with no trauma for the baby. After a caesarean, in addition to «skin to skin» contact, it takes the hypothermia trauma out.



EASY TO USE:

- Ultra compact design with a sleeping cocoon of 35 x 60cm.
- LCD screen with real time information displayed and « turn & push » knob for improving settings and navigation.
- Proclive position : 0 to 25%.
- Automatic cutaneous sensor management.
- Removable Interior compartment for an easier cleaning and disinfection.
- Convenient cables arrangement under the cot.



COMFORTABLE AND SAFE FOR THE BABY:

- PU sleeping cocoon.
- Twins cocoon also available (COCONTWIN).
- Mellow gel mattress.
- The 8 temperature sensors manage the heating with a priority on the hottest sensor.
- Visual and audible alarms for : temperature limits, running errors, accidental power shut down...



MADE IN FRANCE:

ECM develops and produces patient warming systems since 1991 in its production site of Saint Leonard de Noblat, near Limoges - FRANCE !

DESCRIPTION	REFERENCE
All-in-one compact warming cot BABYLEO VITA 48	BABYLEO VITA48

Accessoires

DESCRIPTION	REFERENCE
Replacement cocoon for BABYLEO VITA 48 cot	COCON
Twin cocoon for BABYLEO VITA 48 cot	COCONTWIN

ECM FRENCH QUALITY

CE MARKING 0459

Warming cots are class IIb medical devices. They comply with current European regulations.

AFTER-SALES SERVICE

ECM provides maintenance and repair services for warming cots.
Contact our after-sales service on : +33 5 55 56 32 40.

ECM WARRANTY

ECM warming cots are guaranteed for 12 months (under normal conditions of use). The warranty does not cover maintenance faults made by the user.



ECM - Electro Concept,
A Spengler Group Company
4 Bd. Henri Barbusse, 87400, Saint Léonard de Noblat, FRANCE

Tel. : +33 5 55 56 32 40
Mail : contact@electroconceptmedical.fr
Web : www.electroconceptmedical.fr