

ECM, expert in hypothermia prevention and treatment !



FABRIQUÉ EN
FRANCE
MADE IN



- Easy to use
- Safe
- Made In France

Clinical studies demonstrated that perioperative hypothermia may cause complications for patients. CB48 heating blankets are patient active warming systems dedicated to operative rooms.



CB48 Heating Blankets for Operating Room

CB48 blankets are not using pulsed air technology but an integrated circuit in a cartridge (prevention of noises and particles). Costs are controlled as they are reusable and do not need consumables. Blankets are controlled by a BCB48 electronic control unit..



EASY TO USE:

- Several sizes available with the same control unit box (BCB48).

	<i>Width</i>	<i>Length</i>
CB48-2 : Pediatric	80cm	100cm
CB48-3 : Upper body	50cm	150cm
CB48-11 : Full body	100cm	150cm

- Select temperature from 30 to 40°C with the BCB48 knob and switch on the machine.
- The blankets are made with PU cover, which is very easy to clean and disinfect after use.



SAFE:

- Precision and safety with 8 independant temperature sensors.
- Equipped with visual and audible alarms to call nurses in case of:
 - Heating element or blanket conductor breaking
 - Temperature sensor default
 - Temperature variation too high
 - Temperature above 43°C
 - Power cutage
 - Main board default
- Low voltage direct current power: 48V, to prevent from monitors interferences.



MADE IN FRANCE:

ECM develops and produces patient warming systems since 1991 in its production site of Saint Leonard de Noblat, near Limoges - FRANCE !





CB48 Heating Blankets for Operating Room

References

DESCRIPTION	REFERENCE
OR Heating Blanket - Pediatric	CB48-2+BCB48
OR Heating Blanket - Adult (half body)	CB48-3+BCB48
OR Heating Blanket - Adult (full body)	CB48-11+BCB48

Control unit box BCB48 and blankets are also available separately



Accessories:

DESCRIPTION	REFERENCE
Cover holder for serum stand	SUPPS



ECM FRENCH QUALITY

CE MARKING 0459

Warming cots are class IIb medical devices. They comply with current European regulations.

AFTER-SALES SERVICE

ECM provides maintenance and repair services for heating blankets.
Contact our after-sales service on : +33 5 55 56 32 40.

ECM WARRANTY

ECM heating blankets are guaranteed for 12 months (under normal conditions of use). The warranty does not cover maintenance faults made by the user.



ECM - Electro Concept,
A Spengler Group Company
4 Bd. Henri Barbusse, 87400, Saint Léonard de Noblat, FRANCE

Tel. : +33 5 55 56 32 40
Mail : contact@electroconceptmedical.fr
Web : www.electroconceptmedical.fr