

ECM, expert in hypothermia prevention and treatment !



- Easy to use
- Safe
- Made In France

Clinical studies demonstrated that perioperative hypothermia may cause complications for patients. CR heating blankets are patient active warming systems dedicated to recovery rooms after surgery.



CR Heating blankets for recovery room

CR heated blankets don't use forced-air technology, but have an integrated circuit in a cartridge (prevention of noise and particles). Costs are controlled as they are reusable and don't need consumables. They aim at preventing the risk of hypothermia for patients who just had surgery, and they are used in hospital rooms for post surgical care.

The blankets are managed through the electronical control unit box with LCD display and keyboard.



EASY TO USE:

- Several sizes available with the same control unit box (BCR).

Width x Length

CR11 Adult : 100 x 150 cm

CR2 Child: 80 x 100 cm

CR1 Newborn: 40 x 45 cm → warming the baby without breaking the bond with the mother in the maternity ward

- Set the temperature on BCR box LCD display and press OK.
- PU-covered blankets are easy to clean after use, using a surface cleaning and disinfectant solution.



SAFE:

- Precision and safety with 4 independant temperature sensors.
- Equipped with visual and audible alarms to call nurses in case of:
 - Heating element or blanket conductor breaking
 - Temperature sensor default
 - Temperature variation too high
 - Temperature above 43°C
 - Power cutage
 - Main board default



MADE IN FRANCE:

ECM develops and produces patient warming systems since 1991 in its production site of Saint Leonard de Noblat, near Limoges - FRANCE !

DESCRIPTION	REFERENCE
Newborn heating blanket for recovery room	CR1+BCR
Child heating blanket for recovery room	CR2+BCR
Adult heating blanket for recovery room	CR11+BCR

Control unit box BCR and blankets are also available separately

ECM FRENCH QUALITY

CE MARKING 0459

Warming cots are class IIb medical devices. They comply with current European regulations.

AFTER-SALES SERVICE

ECM provides maintenance and repair services for heating blankets. Contact our after-sales service on : +33 5 55 56 32 40.

ECM WARRANTY

ECM heating blankets are guaranteed for 24 months (under normal conditions of use). The warranty does not cover maintenance faults made by the user.



ECM - Electro Concept,
A Spengler Group Company
4 Bd. Henri Barbusse, 87400, Saint Léonard de Noblat, FRANCE

Tel. : +33 5 55 56 32 40
Mail : contact@electroconceptmedical.fr
Web : www.electroconceptmedical.fr