

ECM highly takes care of very little patients, offering its **warming mattress PrimaLéo** dedicated to warm up newborns in **neonatology and maternity department**.



- Comfortable for the baby
- Easy to use and safe
- Made in France



PrimaLéo Warming Mattress

The **PrimaLéo** removable solution was developed to be used with most types of cots. It is composed of a control unit box, the cartridge heating element and the cocoon/nest.



COMFORTABLE FOR THE BABY:

- The cartridge heating element is made of a gel mattress with protection foam for a cosy feeling.
- The very comfortable PU cocoon allows the baby to be perfectly wedged for more comfort and safety.



EASY TO USE AND SAFE:

- Start up and setting of the warming system operate in a few seconds.
- The electronic system of 8 temperature sensors offer heating security and efficiency with only a 40W power (no interference with other monitors)
- A skin temperature sensor included, carries alarms with high and low thresholds. The control unit box screen displays then, the skin temperature of the baby and the one of the mattress.
- Cleaning and disinfection (quaternary ammonium) are facilitated thanks to the waterproof welding of both upper and lower foils of the cocoon.



MADE IN FRANCE:

ECM develops and produces patient warming systems since 1991 in its production site of Saint Leonard de Noblat, near Limoges - FRANCE !

DESCRIPTION	REFERENCE
Warming mattress	PRIMALEO

ECM FRENCH QUALITY

CE MARKING 0459

Warming mattresses are class IIb medical devices. They comply with current European regulations.

AFTER-SALES SERVICE

ECM provides maintenance and repair services for heated mattresses. Contact our after-sales service on : +33 5 55 56 32 40.

ECM WARRANTY

ECM Warming mattresses are guaranteed for 12 months (under normal conditions of use). The warranty does not cover maintenance faults made by the user.



ECM - Electro Concept,
A Spengler Group Company
4 Bd. Henri Barbusse, 87400, Saint Léonard de Noblat, FRANCE

Tel. : +33 5 55 56 32 40
Mail : contact@electroconceptmedical.fr
Web : www.electroconceptmedical.fr