

The **BABY LEO VITA 48** allows all new borns to get the ideal temperature and warms up the « low weight » ones non able to maintain their necessary normothermia the first days of their life. BABY LEO VITA 48 warming cot is dedicated to neonatology and maternity departments.



- Easy to use
- Comfortable and safe for the baby
- Made in France !



BABY LEO VITA 48 Warming Cot

TECHNICAL SPECIFICATIONS

WEIGHT	17 Kg	
MATERIAL	Mattress	Cross-linked PU gel
	Cocoon	PU
TEMPERATURE CONTROL	Measurement principle	CTN
	Measurement range	20°C to 41°C
	Measurement accuracy	+/- 0.8°C
	Adjustment range	30°C à 38.5°C / 0.5°C step
ELECTRICAL SPECIFICATIONS	Operating voltage	230 Vac - 50 Hz
	DC Power Supply	48V
	Mains filter	Medical specific
	Power consumption	100W (Heater 80 W)
	Electrical classification	CLASS I
	Leakage current	< 10 µA
	Fuses	2 x 3.15A FI (dim. 5x20)
	Thermal fuses card	80°C +/- 10°C 2.5A min
	Thermal fuses heats	55°C +/- 5°C 2.5A max
	CONDITIONS OF USE	Recommended room temperature
Temperature range		10°C / 40°C
Atmospheric pressure		700 hPa to 1060 hPa
Hygrometry (non condensing)		30% to 90%
STORAGE CONDITIONS	Type of packaging	Origin
	Local	Ventilated
	Storage temperature	-20°C to 50°C
	Atmospheric pressure	900 hPa to 1060 hPa
	Hygrometry (non condensing)	30% to 90%
CE MARKING	CE Marking 0459, notified body : G-MED, MD Class : IIb	

DESCRIPTION	REFERENCE
All-in-one compact warming cot BABYLEO VITA 48	BABYLEO VITA48

Accessoires

DESCRIPTION	REFERENCE
Replacement cocoon for BABYLEO VITA 48 cot	COCON
Twin cocoon for BABYLEO VITA 48 cot	COCONTWIN

ECM FRENCH QUALITY

CE MARKING 0459

Warming cots are class IIb medical devices. They comply with current European regulations.

AFTER-SALES SERVICE

ECM provides maintenance and repair services for warming cots. Contact our after-sales service on : +33 5 55 56 32 40.

ECM WARRANTY

ECM warming cots are guaranteed for 12 months (under normal conditions of use). The warranty does not cover maintenance faults made by the user.



ECM - Electro Concept,
A Spengler Group Company
4 Bd. Henri Barbusse, 87400, Saint Léonard de Noblat, FRANCE

Tel. : +33 5 55 56 32 40
Mail : contact@electroconceptmedical.fr
Web : www.electroconceptmedical.fr