

ECM, expert in hypothermia prevention and treatment !



- Easy to use
- Safe
- Made In France



CR Heating blankets for recovery room

CONTROL BOX TECHNICAL SPECIFICATIONS

WEIGHT	Box	0.9 Kg
TEMPERATURE CONTROL	Measurement principle	CTN
	Measurement range	20°C to 41°C
	Measurement accuracy	+/- 0.8°C
	Adjustment range	30°C to 40°C / 0.5°C step
CONDITIONS OF USE	Recommended room temperature	21°C to 26°C
	Temperature range	10°C / 40°C
	Atmospheric pressure	700 hPa to 1060 hPa
	Hygrometry (non condensing)	30% to 90%
STORAGE CONDITIONS	Type of packaging	Origin
	Local	Ventilated
	Storage temperature	-20°C to 50°C
	Atmospheric pressure	900 hPa to 1060 hPa
CE MARKING	Hygrometry (non condensing)	30% to 90%
	CE Marking 0459, notified body : G-MED, Device Class : IIb	

HEATING BLANKETS TECHNICAL SPECIFICATIONS

CR HEATING BLANKET REFERENCE:

	CR1 Model Newborn	CR2 Model Child	CR11 Model Adult
DIMENSIONS	40 x 45 cm	80 x 100 cm	100 x 150 cm
WEIGHT (w/out box) +/-0.2 Kg	0.6 Kg	1.5 Kg	2.4 Kg
POWER	25 W	65 W	150 W
POWER SUPPLY	230 V 50 Hz		
BOX WATERPROOFING	IP 32		
BLANKET WATERPROOFING	IP 35		
COVER	PU		
INNER FABRIC	83% acrylic / 17% cotton		
HANDLE	0	0	4
TEMPERATURE SETTING	30 à 40°C		
TEMPERATURE CONTROL ACCURACY	+/- 1°C		

DESCRIPTION	REFERENCE
Newborn heating blanket for recovery room	CR1+BCR
Children heating blanket for recovery room	CR2+BCR
Adult heating blanket for recovery room	CR11+BCR

Control unit box BCR and blankets are also available separately

ECM FRENCH QUALITY

CE MARKING 0459

Warming cots are class IIb medical devices. They comply with current European regulations.

AFTER-SALES SERVICE

ECM provides maintenance and repair services for heating blankets. Contact our after-sales service on : +33 5 55 56 32 40.

ECM WARRANTY

ECM heating blankets are guaranteed for 24 months (under normal conditions of use). The warranty does not cover maintenance faults made by the user.



ECM - Electro Concept,
A Spengler Group Company
4 Bd. Henri Barbusse, 87400, Saint Léonard de Noblat, FRANCE

Tel. : +33 5 55 56 32 40
Mail : contact@electroconceptmedical.fr
Web : www.electroconceptmedical.fr